

DSM/DTX 系列变档显微镜

SERIES STEP STEREO MICROSCOPES

FEATURES 性能

- ◆ 视场宽广，成像清晰，立体感强
- ◆ 45° 倾斜双目观察，360° 旋转，锁定目镜，左右目镜筒视度调节 ±5dp，出射瞳距调节范围为55-75mm
- ◆ 90° 旋转式物镜，有三种可供选择：1x/2x, 1x/3x, 2x/4x。DTX10,DTX20物镜倍率变换为手轮式转换
- ◆ 产品通过“CE”认证
- ◆ 输入电压90V~240V宽电压
- ◆ 照明：L1. 照射光12V 10W反射型卤素灯，透射光12V 10W卤素灯，亮度可调节。
L2. 照射光和透射光均用LED灯



DSM



DTX10-W



DTX20-W

- ◆ Sharp stereo erect image over a wide field of view.
- ◆ 45° inclined and 360° rotatable binocular head with locked-in eyepiece. Both ocular-tube with diopter adjustment ±5dp. The interpupillary distance is between 55-75mm.
- ◆ Three kind of 90° rotatable objective for selection:1x/2x,1x/3x,2x/4x. DTX10/DTX20 with knob rotation.
- ◆ Product has CE approval.
- ◆ Power supply 90V~240V.
- ◆ Illumination: L1. Incident reflected halogen light 12V10W;
Transmitted halogen light 12V10W,Brigtness is adjustable.
L2. Both incident and transmitted light with LED.



DSMW



DTX11-W



DTX21

SPECIFICATION 规格

| 型号 Model | 目镜 Eyepiece | 物镜 Objective | 升降机构 Up & down Mechanism | 底座/工作台 Base/Working Stage |
|-------------|----------------|-----------------|-----------------------------|--|
| DSM | | | 立柱式 Pillar style | Ø95mm 黑白工作台 Ø95mm White/Black working stage 扇形高平台底座 Fan high base |
| DSM-W | | 1X-2X | 弯臂式 Arm curved style | Ø95mm 黑白工作台 Ø95mm White/Black working stage 扇形高平台底座 Fan high base |
| DTX10-W | SWF10X | 1X-3X | 弯臂式 Arm curved style | Ø95mm 黑白工作台 Ø95mm White/Black working stage 扇形平底座 Fan plain base |
| DTX11-W | | | 弯臂式 Arm curved style | Ø95mm 磨砂玻璃 Ø95mm Frosted glass 扇形平底座 Fan plain base |
| DTX20-W | | 2X-4X | 弯臂式 Arm curved style | Ø95mm 黑白工作台 Ø95mm White/Black working stage 扇形高平台底座 Fan high base |
| DTX21 | | | 立柱式 Pillar style | Ø95mm 磨砂玻璃 Ø95mm Frosted glass 扇形高平台底座 Fan high base |

可选购附件 (Optional accessories): WF5x、SWF10x、WF15x、WF20x 目镜 (Eyepiece); 环形灯 (Ring light); 0.5X、0.75X、1.5X、2X 附加物镜 (Auxiliary Objective)

OPTICAL DATA 光学参数

| 物镜放大率 Objective Magnification | 目镜 Eyepiece | | | | | | | | 工作距离 Working Distance (mm) |
|-------------------------------------|----------------------|-----------------------------|----------------------|-----------------------------|----------------------|-----------------------------|----------------------|-----------------------------|----------------------------------|
| | WF5X | | SWF10X | | WF15X | | WF20X | | |
| | 放大率 Magnification | 视场 Field of View (mm) | 放大率 Magnification | 视场 Field of View (mm) | 放大率 Magnification | 视场 Field of View (mm) | 放大率 Magnification | 视场 Field of View (mm) | |
| 1X | 5X | 20.0 | 10X | 22.0 | 15X | 15.0 | 20X | 10.0 | 100 |
| 2X | 10X | 10.0 | 20X | 11.0 | 30X | 7.5 | 40X | 6.5 | 100 |
| 1X | 5X | 20.0 | 10X | 22.0 | 15X | 15.0 | 20X | 10.0 | 100 |
| 3X | 15X | 7.0 | 15X | 7.3 | 45X | 5.0 | 60X | 4.3 | 100 |
| 2X | 10X | 10.0 | 20X | 11.0 | 30X | 7.5 | 40X | 6.5 | 100 |
| 4X | 20X | 5.0 | 40X | 5.5 | 60X | 3.7 | 80X | 3.2 | 100 |

ZTX-10/20 系列变档显微镜

SERIES STEP STEREO MICROSCOPES

FEATURES 性能

- ◆ 视场宽广，成像清晰，立体感强
- ◆ 45° 倾斜双目观察，360° 旋转，锁定目镜，左目镜筒视度调节 ± 5dp，出射瞳距调节范围为55-75mm
- ◆ 90° 旋转式物镜，有三种可供选择：1x/2x, 1x/3x, 2x/4x
- ◆ ZTX-10和ZTX-20系列为立柱式支架，ZTX-10-C型为扇形底座，ZTX-20-W系列为弯臂式支架
- ◆ 产品通过CE认证
- ◆ 输入电压为110V/60Hz, 220V/50Hz, 90V~240V宽电压
- ◆ 照明：L1. 照射光12V 10W反射型卤素灯，透射光12V 10W卤素灯
L2. 照射光12V 10W反射型卤素灯，透射光12V 10W卤素灯，亮度可调节
L3. 照射照明12V 10W卤素灯，透射照明5W荧光灯
L4. 照射照明12V 10W卤素灯，透射照明5W荧光灯，亮度可调节
L5. 透射及照射照明均用LED灯
- ◆ 多种附件可供选择
- ◆ 每台泡沫包装，放入内纸箱，四台一个外箱，外包装尺寸为ZTX-10/20: 85x35x40cm ZTX-20-W: 50x37x48cm



ZTX-10



ZTX-20

- ◆ Sharp stereo erect image over a wide field of view.
- ◆ 45° inclined and 360° rotatable binocular head with locked-in eyepiece. The left ocular-tube with diopter adjustment ±5dp. The interpupillary distance is between 55-75mm.
- ◆ Three kind of 90° rotatable objective for selection: 1x/2x, 1x/3x, 2x/4x
- ◆ ZTX-10/ZTX-20 series supplied with pillar style stand; ZTX-10-C series supplied with fan base and pillar stand, ZTX-20-W series supplied with arm curved style stand.
- ◆ Product has CE approval.
- ◆ Power supply 110V/60Hz or 220V/50Hz or 90V~240V.
- ◆ Illumination: L1. Incident reflected halogen lamp 12V10W; Transmitted halogen lamp 12V10W.
L2. Incident reflected halogen lamp 12V10W; Transmitted halogen lamp 12V10W; Brightness is adjustable.
L3. Halogen lamp 12V/10W for incident lamp and fluorescent lamp 5W for transmitted lamp.
L4. Halogen lamp 12V/10W for incident lamp and fluorescent lamp 5W for transmitted lamp. Brightness is adjustable.
L5. Both incident and transmitted light with LED.
- ◆ A wide range of accessories for selection.
- ◆ Each set packed in styrofoam and in an inner box. Four sets in one carton. The outer measurement is ZTX-10/20: 85x35x40cm ZTX-20-W: 50x37x48cm.

SPECIFICATION 规格

| 型号 Model | 目镜 Eyepiece | 物镜 Objective | 工作台 Working Stage | 照明 Illumination |
|-----------------------|-------------|--------------|--------------------|--|
| ZTX-10 | WF10X | | 黑/白 White/black | 自然光 Natural illumination |
| ZTX-10-W | WF10X | Rotated | 黑/白 White/black | 自然光 Natural illumination |
| ZTX-10-C ₁ | WF10X | 1x-2x | 黑/白 White/black | 自然光 Natural illumination |
| ZTX-20 | WF10X | 1x-3x | 磨砂玻璃 Frosted glass | 照射光/透射光 Incident illumination & transmitted illumination |
| ZTX-20-W | WF10X | 2x-4x | 磨砂玻璃 Frosted glass | 照射光/透射光 Incident illumination & transmitted illumination |
| ZTX-20-C ₂ | WF10X | | 磨砂玻璃 Frosted glass | 照射光/透射光 Incident illumination & transmitted illumination |

可选购附件 (Optional accessories): 0.5x、0.75x、1.5x、2x 附加物镜 (Auxiliary Objective); WF5x、WF10x、WF15x、WF20x 目镜 (Eyepiece); 环形灯 (Ring light)



ZTX-20-C₂



ZTX-20-E

OPTICAL DATA 光学参数

| 目镜 Eyepiece | 标准物镜 Standard Objective | | | 0.5X | | 附加物镜 Auxiliary Objective 0.75X | | 1.5X | | 2.0X | |
|-------------|-------------------------|-------------------------|----------------------------|-------------------------|----------------------------|--------------------------------|----------------------------|-------------------------|----------------------------|-------------------------|----------------------------|
| | 物镜倍率 Objective | 视场范围 Field of View (mm) | 工作距离 Working Distance (mm) | 视场范围 Field of View (mm) | 工作距离 Working Distance (mm) | 视场范围 Field of View (mm) | 工作距离 Working Distance (mm) | 视场范围 Field of View (mm) | 工作距离 Working Distance (mm) | 视场范围 Field of View (mm) | 工作距离 Working Distance (mm) |
| WF5X | 1X | 20 | 109 | 37.5 | 120 | 24 | 75 | 14 | 40 | 11 | 25 |
| | 2X | 10 | 99 | 19 | 120 | 12 | 75 | 7 | 40 | 5.5 | 25 |
| | 3X | 7 | 84 | 12.5 | 120 | 8 | 75 | 4.5 | 40 | 3.5 | 25 |
| | 4X | 5 | 81 | 19 | 120 | 6 | 75 | 3.5 | 40 | 2.8 | 25 |
| WF10X | 1X | 20 | 109 | 37.5 | 120 | 25 | 75 | 14 | 40 | 11 | 25 |
| | 2X | 10 | 99 | 19 | 120 | 12 | 75 | 7 | 40 | 5.5 | 25 |
| | 3X | 7 | 84 | 12.5 | 120 | 8 | 75 | 4.5 | 40 | 3.5 | 25 |
| WF15X | 1X | 15 | 109 | 30 | 120 | 18 | 75 | 14 | 40 | 8 | 25 |
| | 2X | 7.5 | 99 | 15 | 120 | 9 | 75 | 7 | 40 | 4 | 25 |
| | 3X | 5 | 84 | 9.5 | 120 | 6 | 75 | 4.5 | 40 | 3 | 25 |
| | 4X | 3.8 | 81 | 7 | 120 | 4.5 | 75 | 3.5 | 40 | 2 | 25 |
| WF20X | 1X | 10 | 109 | 19 | 120 | 12 | 75 | 7 | 40 | 5.5 | 25 |
| | 2X | 5 | 99 | 10 | 120 | 6 | 75 | 3.5 | 40 | 2.5 | 25 |
| | 3X | 3.5 | 84 | 6.2 | 120 | 4 | 75 | 2.4 | 40 | 2 | 25 |
| | 4X | 2.5 | 81 | 4.7 | 120 | 3 | 75 | 2 | 40 | 1.5 | 25 |

XTX-1/4 系列变档显微镜

SERIES STEP STEREO MICROSCOPES

FEATURES 性能

- ◆ 视场宽广，成像清晰，立体感强
- ◆ 垂直双目观察，锁定目镜，左目镜筒视度调节±5dp，出射瞳距调节范围为55-75mm
- ◆ XTX-1 V型槽插入式物镜、带锁定机构、便于更换；物镜1x, 2x, 3x, 4x, 6x可选择
- ◆ XTX-4型固定式物镜2X
- ◆ 产品通过CE认证
- ◆ 输入电压为110V/60Hz或220V/50Hz
- ◆ 目镜：WF5x WF10x WF15x WF20x可选
- ◆ 照明：12V 10W 白炽灯
- ◆ 每台泡沫包装，放入内纸箱，四台一个外箱，外包装尺寸：68x31x37cm



XTX-1A



XTX-1B



XTX-1C

SPECIFICATION 规格

| 型号 Model | 目镜 Eyepiece | 物镜 Objective | 工作台 Working Stage | 照明 Illumination |
|----------|-------------|--------------|--------------------|--|
| XTX-1A | WF10X | Inserted 2X | 黑/白 White/black | 自然光 Natural illumination |
| XTX-1B | WF10X | Inserted 2X | 黑/白 White/black | 照射光 Incident illumination |
| XTX-1C | WF10X | Inserted 2X | 磨砂玻璃 Frosted glass | 照射光/透射光 Incident illumination & transmitted illumination |
| XTX-4A | WF10X | Fixed 2X | 黑/白 White/black | 自然光 Natural illumination |
| XTX-4B | WF10X | Fixed 2X | 黑/白 White/black | 照射光 Incident illumination |
| XTX-4C | WF10X | Fixed 2X | 磨砂玻璃 Frosted glass | 照射光/透射光 Incident illumination & transmitted illumination |



XTX-4A



XTX-4B



XTX-4C

- ◆ Sharp stereo erect image over a wide field of view.
- ◆ Upright binocular head with locked-in eyepiece. The left ocular-tube with diopter adjustment±5dp. The interpupillary distance is between 55-75mm.
- ◆ XTX-1 series V-shape inserted style objective with locked-in device ensure easy replacement. Objective: 1x, 2x, 3x, 4x, 6x for selection.
- ◆ XTX-4 series with fixed style objective 2X.
- ◆ Product has CE approval.
- ◆ Power supply 110V/60Hz or 220V/50Hz.
- ◆ Eyepiece: WF5x WF10x WF15x WF20x for selection.
- ◆ Illumination: 12V 10W incandescent lamp.
- ◆ Each set packed in styrofoam and in an inner box. Four sets in one carton.The outer measurement is 68x31x37cm.

OPTICAL DATA 光学参数

| 物镜放大率 Objective Magnification | 目镜 Eyepiece | | | | | | | | 工作距离 Working Distance (mm) |
|-------------------------------|-------------------|-----------------------|-------------------|-----------------------|-------------------|-----------------------|-------------------|-----------------------|----------------------------|
| | WF5X | | SWF10X | | WF15X | | WF20X | | |
| | 放大率 Magnification | 视场 Field of View (mm) | 放大率 Magnification | 视场 Field of View (mm) | 放大率 Magnification | 视场 Field of View (mm) | 放大率 Magnification | 视场 Field of View (mm) | |
| 1X | 5X | 20.0 | 10X | 20.0 | 15X | 15.0 | 20X | 10.0 | 57 |
| 2X | 10X | 10.0 | 20X | 10.0 | 30X | 7.5 | 40X | 6.5 | 80 |
| 3X | 15X | 7.0 | 30X | 6.7 | 45X | 5.0 | 60X | 4.3 | 61 |
| 4X | 20X | 5.0 | 40X | 5.0 | 60X | 3.7 | 80X | 3.2 | 57 |
| 6X | 30X | 4.0 | 60X | 3.3 | 90X | 2.5 | 120X | 2.1 | 66 |

XTX-2/3 系列变档显微镜

SERIES STEP STEREO MICROSCOPES

FEATURES 性能

- ◆ 视场宽广，成像清晰，立体感强
- ◆ 45° 倾斜双目观察，锁定目镜，左目镜筒视度调节 ± 5dp，出射瞳距调节范围为55-75mm
- ◆ XTX-2系列V型槽插入式物镜、带锁定机构、便于更换，物镜1x, 2x, 3x, 4x, 6x可选择
- ◆ XTX-3系列180° 旋转式物镜，有三种可供选择：1x/2x, 1x/3x, 2x/4x
- ◆ 产品通过CE认证
- ◆ 输入电压为110V/60Hz或220V/50Hz
- ◆ 目镜：WF5x WF10x WF15x WF20x可选
- ◆ 照明：1. 白炽灯12V 10W 2. 卤素灯12V 10W 3. 高亮度LED灯
- ◆ 每台泡沫包装，放入内纸箱，四台一个外箱。外包装尺寸：68x31x37cm



XTX-2B



XTX-2C



XTX-2A

- ◆ Sharp stereo erect image over a wide field of view.
- ◆ 45° inclined head with locked-in eyepiece. The left ocular-tube with diopter adjustment ± 5dp. The interpupillary distance is between 55-75mm.
- ◆ XTX-2 series V-shape inserted style objective with locked-in device ensure easy replacement. Objective: 1x, 2x, 3x, 4x, 6x for selection.
- ◆ XTX-3 series has three kind of 180° rotatable objective for selection: 1x/2x, 1x/3x, 2x/4x.
- ◆ Product has CE approval.
- ◆ Power supply 110V/60Hz or 220V/50Hz.
- ◆ Eyepiece: WF5x WF10x WF15x WF20x for selection.
- ◆ Illumination: 1. Incandescent lamp 12V 10W 2. Halogen lamp 12V 10W 3. High brightness LED lamp.
- ◆ Each set packed in styrofoam and in an inner box. Four sets in one carton. The outer measurement is 68x31x37cm.

SPECIFICATION 规格

| 型号 Model | 目镜 Eyepiece | 物镜 Objective | 工作台 Working Stage | 照明 Illumination |
|----------|-------------|--------------|--------------------|--|
| XTX-2A | WF10X | Inserted 2X | 黑/白 White/black | 自然光 Natural illumination |
| XTX-2B | WF10X | Inserted 2X | 黑/白 White/black | 照射光 Incident illumination |
| XTX-2C | WF10X | Inserted 2X | 磨砂玻璃 Frosted glass | 照射光 / 透射光 Incident illumination & transmitted illumination |
| XTX-3A | WF10X | Rotated | 黑/白 White/black | 自然光 Natural illumination |
| XTX-3A-W | WF10X | 1X-2X | 黑/白 White/black | 自然光 Natural illumination |
| XTX-3B | WF10X | 1X-3X | 黑/白 White/black | 照射光 Incident illumination |
| XTX-3C | WF10X | 2X-4X | 磨砂玻璃 Frosted glass | 照射光 / 透射光 Incident illumination & transmitted illumination |
| XTX-3C-W | WF10X | 2X-4X | 磨砂玻璃 Frosted glass | 照射光 / 透射光 Incident illumination & transmitted illumination |



XTX-3A



XTX-3B



XTX-3C

OPTICAL DATA 光学参数

| 型号 Model | 物镜放大率 Objective Magnification | 目镜 Eyepiece | | | | | | | | 工作距离 Working Distance (mm) |
|----------|-------------------------------|-----------------------|-------------------|-----------------------|-------------------|-----------------------|-------------------|-----------------------|------|----------------------------|
| | | WF5X | | SWF10X | | WF15X | | WF20X | | |
| | 放大率 Magnification | 视场 Field of View (mm) | 放大率 Magnification | 视场 Field of View (mm) | 放大率 Magnification | 视场 Field of View (mm) | 放大率 Magnification | 视场 Field of View (mm) | | |
| XTX-2A | 2X | 10X | 10.0 | 20X | 10.0 | 30X | 7.5 | 40X | 6.5 | 80 |
| XTX-2B | 3X | 15X | 7.0 | 30X | 6.7 | 45X | 5.0 | 60X | 4.3 | 61 |
| XTX-2C | 4X | 20X | 5.0 | 40X | 5.0 | 60X | 3.7 | 80X | 3.2 | 57 |
| | 6X | 30X | 4.0 | 60X | 3.3 | 90X | 2.5 | 120X | 2.1 | 66 |
| XTX-3A | 1X | 5X | 20.0 | 10X | 20.0 | 15X | 15.0 | 20X | 10.0 | 57 |
| | 2X | 10X | 10.0 | 20X | 10.0 | 30X | 7.5 | 40X | 6.5 | 57 |
| XTX-3B | 1X | 5X | 20.0 | 10X | 20.0 | 15X | 15.0 | 20X | 10.0 | 57 |
| | 3X | 15X | 7.0 | 30X | 6.7 | 45X | 5.0 | 60X | 4.3 | 57 |
| XTX-3C | 2X | 10X | 10.0 | 20X | 10.0 | 30X | 7.5 | 40X | 6.5 | 57 |
| | 4X | 20X | 5.0 | 40X | 5.0 | 60X | 3.7 | 80X | 3.2 | 57 |

XTX-5 系列变档显微镜

SERIES STEP STEREO MICROSCOPES

FEATURES 性能

- ◆ 视场宽广，成像清晰，立体感强
- ◆ 45° 倾斜双目观察，锁定目镜，头部可360° 转动，左目镜筒视度调节 ± 5dp，出射瞳距调节范围为55-75mm
- ◆ 90° 旋转式物镜，有三种可供选择：1x/2x, 1x/3x, 2x/4x
- ◆ XTX-5A/5B/5C系列为立柱式支架，XTX-5A-W,XTX-5C-W系列为弯臂式支架
- ◆ 产品通过CE认证
- ◆ 输入电压为110V/60Hz或220/50Hz
- ◆ 照明：1.白炽灯12V 10W 2. 卤素灯12V 10W 3. 高亮度LED灯
- ◆ 每台泡沫包装，放入内纸箱，四台一个外箱，外包装尺寸为：68x31x37cm



XTX-5A



XTX-5B



XTX-5C

SPECIFICATION 规格

| 型号 Model | 目镜 Eyepiece | 物镜 Objective | 工作台 Working Stage | 照明 Illumination |
|----------|-------------|--------------|--------------------|--|
| XTX-5A | WF10X | Rotated | 黑 / 白 White/black | 自然光 Natural illumination |
| XTX-5B | WF10X | 1X-2X | 黑 / 白 White/black | 照射光 Incident illumination |
| XTX-5C | WF10X | 1X-3X | 磨砂玻璃 Frosted glass | 照射光 / 透射光 Incident illumination & transmitted illumination |
| XTX-5A-W | WF10X | 2X-4X | 黑 / 白 White/black | 自然光 Natural illumination |
| XTX-5C-W | WF10X | | 磨砂玻璃 Frosted glass | 照射光 / 透射光 Incident illumination & transmitted illumination |



XTX-5A-W



XTX-5C-W

OPTICAL DATA 光学参数

| 物镜放大率 Objective Magnification | 目镜 Eyepiece | | | | | | | | 工作距离 Working Distance (mm) |
|-------------------------------|-------------------|-----------------------|-------------------|-----------------------|-------------------|-----------------------|-------------------|-----------------------|----------------------------|
| | WF5X | | SWF10X | | WF15X | | WF20X | | |
| | 放大率 Magnification | 视场 Field of View (mm) | 放大率 Magnification | 视场 Field of View (mm) | 放大率 Magnification | 视场 Field of View (mm) | 放大率 Magnification | 视场 Field of View (mm) | |
| 1X | 5X | 20.0 | 10X | 20.0 | 15X | 15.0 | 20X | 10.0 | 76 |
| 2X | 10X | 10.0 | 20X | 10.0 | 30X | 7.5 | 40X | 6.5 | 76 |
| 1X | 5X | 20.0 | 10X | 20.0 | 15X | 15.0 | 20X | 10.0 | 76 |
| 3X | 15X | 7.0 | 30X | 6.7 | 45X | 5.0 | 60X | 4.3 | 76 |
| 2X | 10X | 10.0 | 20X | 10.0 | 30X | 7.5 | 40X | 6.5 | 76 |
| 4X | 20X | 5.0 | 40X | 5.0 | 60X | 3.7 | 80X | 3.2 | 76 |

- ◆ Sharp stereo erect image over a wide field of view.
- ◆ 45° inclined and 360° rotatable binocular head with locked-in eyepiece. The left ocular-tube with diopter adjustment ±5dp. The interpupillary distance is between 55-75mm.
- ◆ Three kind of 90° rotatable objective for selection: 1x/2x, 1x/3x, 2x/4x.
- ◆ XTX-5A/5B/5C series supplied with pillar style stand: XTX-5A-W,XTX-5C-W series supplied with arm curved style stand.
- ◆ Product has CE approval.
- ◆ Power supply 110V/60Hz or 220V/50Hz.
- ◆ Illumination: 1.Incandescent lamp12V 10W 2.Halogen lamp12V 10W 3. High brightness LED lamp.
- ◆ Each set packed in styrofoam and in an inner box. Four sets in one carton.The outer measurement is 68x31x37cm.

XTX-6/7 系列变档显微镜

SERIES STEP STEREO MICROSCOPES

FEATURES 性能

- ◆ 视场宽广，成像清晰，立体感强
- ◆ 45° 倾斜双目观察，锁定目镜，头部可360° 转动，左目镜筒视度调节 ±5dp，出射瞳距调节范围为55-75mm
- ◆ 90° 旋转式物镜，有三种倍率可供选择：1x/2x, 1x/3x, 2x/4x
- ◆ XTX-6A,XTX-6C系列为立柱式支架，XTX-6A-W,XTX-6C-W,XTX-7A-W,XTX-7C-W系列为弯臂式支架
- ◆ 产品通过CE认证
- ◆ 输入电压为110V/60Hz或220/50Hz
- ◆ 照明：1. 白炽灯12V 10W 2. 卤素灯12V 10W 3. 高亮度LED灯
- ◆ 每台泡沫包装，放入内纸箱，四台一个外箱，外包装尺寸为：68x31x37cm



XTX-6A



XTX-6C



XTX-6C-W

SPECIFICATION 规格

| 型号 Model | 目镜 Eyepiece | 物镜 Objective | 工作台 Working Stage | 照明 Illumination |
|----------|-------------|--------------|--------------------|--|
| XTX-6A | WF10X | | 黑/白 White/black | 自然光 Natural illumination |
| XTX-6C | WF10X | Rotated | 磨砂玻璃 Frosted glass | 照射光/透射光 Incident illumination & transmitted illumination |
| XTX-6A-W | WF10X | 1X-2X | 黑/白 White/black | 自然光 Natural illumination |
| XTX-6C-W | WF10X | 1X-3X | 磨砂玻璃 Frosted glass | 照射光/透射光 Incident illumination & transmitted illumination |
| XTX-7A-W | WF10X | 2X-4X | 磨砂玻璃 Frosted glass | 自然光 Natural illumination |
| XTX-7C-W | WF10X | | 磨砂玻璃 Frosted glass | 照射光/透射光 Incident illumination & transmitted illumination |



XTX-6A-W



XTX-7A-W



XTX-7C-W

OPTICAL DATA 光学参数

| 物镜放大率 Objective Magnification | 目镜 Eyepiece | | | | | | | | 工作距离 Working Distance (mm) |
|-------------------------------|-------------------|-----------------------|-------------------|-----------------------|-------------------|-----------------------|-------------------|-----------------------|----------------------------|
| | WF5X | | SWF10X | | WF15X | | WF20X | | |
| | 放大率 Magnification | 视场 Field of View (mm) | 放大率 Magnification | 视场 Field of View (mm) | 放大率 Magnification | 视场 Field of View (mm) | 放大率 Magnification | 视场 Field of View (mm) | |
| 1X | 5X | 20.0 | 10X | 20.0 | 15X | 15.0 | 20X | 10.0 | 76 |
| 2X | 10X | 10.0 | 20X | 10.0 | 30X | 7.5 | 40X | 6.5 | 76 |
| 1X | 5X | 20.0 | 10X | 20.0 | 15X | 15.0 | 20X | 10.0 | 76 |
| 3X | 15X | 7.0 | 30X | 6.7 | 45X | 5.0 | 60X | 4.3 | 76 |
| 2X | 10X | 10.0 | 20X | 10.0 | 30X | 7.5 | 40X | 6.5 | 76 |
| 4X | 20X | 5.0 | 40X | 5.0 | 60X | 3.7 | 80X | 3.2 | 76 |

- ◆ Sharp stereo erect image over a wide field of view.
- ◆ 45° inclined and 360° rotatable binocular head with locked-in eyepiece. The left ocular-tube with diopter adjustment ±5dp. The interpupillary distance is between 55-75mm.
- ◆ Three kinds of 90° rotational objective for selection: 1x/2x, 1x/3x, 2x/4x.
- ◆ XTX-6A,XTX-6C series supplied with pillar style stand: XTX-6A-W,XTX-6C-W,XTX-7A-W,XTX-7C-W series supplied with arm curved style stand.
- ◆ Product has CE approval.
- ◆ Power supply 110V/60Hz or 220V/50Hz.
- ◆ Illumination: 1.Incandescent lamp12V 10W 2.Halogen lamp12V 10W 3. High brightness LED lamp.
- ◆ Each set packed in styrofoam and in an inner box. Four sets in one carton.The outer measurement is 68x31x37cm.

XTX-8/9 系列变档显微镜

SERIES STEP STEREO MICROSCOPES



XTX-8C-2



XTX-8C-3

SPECIFICATION 规格

| 型号 Model | 目镜 Eyepiece | 物镜 Objective | 工作台 Working Stage | 照明 Illumination |
|-------------|----------------|-------------------------|-----------------------------------|---|
| XTX-8C-2 | WF10X | 1X-2X 1X-3X 2X-4X | Ø75mm 磨砂玻璃 Ø75mm Frosted glass | 照射光 / 透射光 Incident illumination & transmitted illumination |
| XTX-8C-3 | WF10X | | Ø75mm 磨砂玻璃 Ø75mm Frosted glass | 照射光 / 透射光 Incident illumination & transmitted illumination |
| XTX-9A | WF10X | 1X-2X-3X 1X-2X-4X | Ø95mm 黑 / 白板 Ø95mm White/black | 照射光 / 透射光 Incident illumination & transmitted illumination |
| XTX-9C | WF10X | | Ø95mm 磨砂玻璃 Ø95mm Frosted glass | 照射光 / 透射光 Incident illumination & transmitted illumination |



XTX-9A



XTX-9C

FEATURES 性能

- ◆ 视场宽广，成像清晰，立体感强
- ◆ 45° 倾斜双目观察，锁定目镜，左目镜筒视度调节 ± 5dp，出射瞳距调节范围为55-75mm
- ◆ 旋转式三组物镜1x/2x/3x或1x/2x/4x
- ◆ XTX-8C系列为弯臂式支架，XTX-9A,XTX-9C为立柱式支架
- ◆ 产品通过CE认证
- ◆ 输入电压为110V/60Hz或220V/50Hz
- ◆ 照明：1. 白炽灯12V 10W 2. 卤素灯12V 10W 3. 高亮度LED灯
- ◆ 每台泡沫包装，放入内纸箱，四台一个外箱，外包装尺寸为：68x31x37cm

- ◆ Sharp stereo erect image over a wide field of view.
- ◆ 45° inclined binocular head with locked-in eyepiece. The left ocular-tube with diopter adjustment ±5dp. The interpupillary distance is between 55-75mm.
- ◆ Three kinds of rotational objective for selection: 1x/2x/3x or 1x/2x/4x.
- ◆ XTX-8C series supplied with arm curved style stand: XTX-9A, XTX-9C supplied with pillar style stand.
- ◆ Product has CE approval.
- ◆ Power supply 110V/60Hz or 220V/50Hz.
- ◆ Illumination: 1.Incandescent lamp12V 10W 2.Halogen lamp12V 10W 3. High brightness LED lamp.
- ◆ Each set packed in styrofoam and in an inner box. Four sets in one carton.The outer measurement is 68x31x37cm.

OPTICAL DATA 光学参数

| 型号 Model | 物镜放大率 Objective Magnification | 目镜 Eyepiece | | | | | | | | 工作距离 Working Distance (mm) |
|-------------|----------------------------------|--------------------------|----------------------|--------------------------|----------------------|--------------------------|----------------------|--------------------------|------|-------------------------------|
| | | WF5X | | SWF10X | | WF15X | | WF20X | | |
| | 放大率 Magnification | 视场 Field of View (mm) | 放大率 Magnification | 视场 Field of View (mm) | 放大率 Magnification | 视场 Field of View (mm) | 放大率 Magnification | 视场 Field of View (mm) | | |
| XTX-8C | 1X | 5X | 20.0 | 10X | 20.0 | 15X | 15.0 | 20X | 10.0 | 100 |
| | 2X | 10X | 10.0 | 20X | 10.0 | 30X | 7.5 | 40X | 6.5 | 100 |
| | 3X | 15X | 7.0 | 30X | 6.7 | 45X | 5.0 | 60X | 4.3 | 100 |
| | 4X | 20X | 5.0 | 40X | 5.0 | 60X | 3.7 | 80X | 3.2 | 100 |
| XTX-9 | 1X | 5X | 20.0 | 10X | 20.0 | 15X | 15.0 | 20X | 10.0 | 57 |
| | 2X | 10X | 10.0 | 20X | 10.0 | 30X | 7.5 | 40X | 6.5 | 57 |
| | 3X | 15X | 7.0 | 30X | 6.7 | 45X | 5.0 | 60X | 4.3 | 57 |
| | 4X | 20X | 5.0 | 40X | 5.0 | 60X | 3.7 | 80X | 3.2 | 57 |

XDT/STX/C-D 系列变档显微镜

SERIES STEP STEREO MICROSCOPES

FEATURES 性能

- ◆ 视场宽广，成像清晰，立体感强
- ◆ 45° 倾斜单目观察，锁定目镜
- ◆ 头部可360° 旋转
- ◆ 放大倍率20x
- ◆ 工作距离80mm
- ◆ XDT-1金属结构，XDT-2金属塑料结构。
- ◆ 单台泡沫包装，10台一个大纸箱，外包装尺寸：76x41x26cm

- ◆ Sharp stereo erect image over a wide field of view.
- ◆ 45° inclined monocular head with locked in eyepiece.
- ◆ 360° rotatable head.
- ◆ Total magnification 20x.
- ◆ Working distance 80mm.
- ◆ XDT-1: metal frame; XDT-2: metal-plastic frame.
- ◆ Each set packed in styrofoam and in an inner box. Ten sets in one carton. The outer measurement is 76x41x26cm.



XDT - 1



XDT - 2

FEATURES 性能

- ◆ 视场宽广，成像清晰，立体感强
- ◆ 垂直双目观察，出瞳距调节范围55-75mm，
- ◆ 目镜WF10X,物镜2X
- ◆ 工作距离70mm
- ◆ 内置照射光: 白炽灯或LED灯
- ◆ 电源为二节五号电池
- ◆ 塑料机身，携带方便
- ◆ 单台泡沫包装10台一个大纸箱，外包装尺寸：69.5x36x29cm.

- ◆ Sharp stereo erect image over a wide field of view.
- ◆ Upright binocular head, the interpupillary distance is between 55-75mm.
- ◆ Eyepiece WF10X, Objective 2X.
- ◆ Working distance: 70mm.
- ◆ Built-in incident light: Incandescent lamp or LED.
- ◆ Power 2 units of AA batteries.
- ◆ Plastic body easily to carry.
- ◆ Each set packed in styrofoam and in an inner box. Ten sets in one carton. The outer measurement is 69.5x36x29cm.



STX - 45



STX

FEATURES 性能

- ◆ 视场宽广，成像清晰，立体感强
- ◆ 左目镜筒视度调节±5dp，出射瞳距调节范围55-75mm
- ◆ 物镜: 0.8X
- ◆ 目镜: WF5X, WF10X, WF15X, WF20X可选
- ◆ 工作距离: 145mm
- ◆ 照明: 环形LED灯，荧光环形灯可选

- ◆ Sharp stereo erect image over a wide field of view
- ◆ The left ocular-tube with diopter adjustment ±5dp. The interpupillary distance is between 55-75mm
- ◆ Objective: 0.8x
- ◆ Eyepiece: WF5X, WF10X, WF15X, WF20X for selection
- ◆ Working distance: 145mm
- ◆ Illumination: LED ring light or fluorescent ring light for selection



C - 6D

FEATURES 性能

- ◆ 视场宽广，成像清晰，立体感强
- ◆ 左目镜筒视度调节±5dp，出射瞳距调节范围55-75mm
- ◆ V型槽插入式物镜，带锁定机构，便于更换，物镜 0.5x,1x可选
- ◆ 照明用高亮度LED光源
- ◆ 目镜，WF5X, WF10X, WF15X, WF20X可选

- ◆ Sharp stereo erect image over a wide field of view.
- ◆ The left ocular-tube with diopter adjustment ±5dp. The interpupillary distance is between 55-75mm.
- ◆ V-shape inserted style objective with locked-in device ensures easy replacement, objective, 0.5x,1x for selection.
- ◆ Illumination: High brightness LED lamp.
- ◆ Eyepiece: WF5X, WF10X, WF15X, WF20X for selection.



C - 2D

OPTICAL DATA 光学参数

| 物镜放大率 Objective Magnification | 目 镜 Eyepiece | | | | | | | | 工作距离 Working Distance (mm) |
|-------------------------------------|----------------------|-----------------------------|----------------------|-----------------------------|----------------------|-----------------------------|----------------------|-----------------------------|----------------------------------|
| | WF5X | | SWF10X | | WF15X | | WF20X | | |
| | 放大率 Magnification | 视场 Field of View (mm) | 放大率 Magnification | 视场 Field of View (mm) | 放大率 Magnification | 视场 Field of View (mm) | 放大率 Magnification | 视场 Field of View (mm) | |
| 0.5X | 2.5X | 40.0 | 5X | 40.0 | 7.5X | 30.0 | 10X | 20.0 | 250 |
| 1X | 5X | 20.0 | 10X | 20.0 | 15X | 15.0 | 20X | 10.0 | 230 |